

pydicom tutorial

January 27, 2022

```
[1]: import pydicom
import matplotlib.pyplot as plt
import os
```

We can use the pydicom function `dcmread` to load a dataset as a pydicom-specific class. Printing the dataset will display all its available attributes, corresponding to Dicom header tags. These may be used to access the data. For example, the image itself is saved as the `pixel_array` attribute. It is read in as a 2D numpy array of unsigned int16 values. Most other attributes are read in as strings.


For reference, here is a basic example.

```
[2]: # first, read in the data and get a list of files
imd = 'imgs/'
L = sorted(os.listdir(imd))

# read the first dataset
i = 0
ds = pydicom.dcmread(imd+L[i])

# extract some key Dicom information
M = ds.pixel_array # the image
rows, cols = ds.Rows, ds.Columns # image dimensions
pir = ds.PixelIntensityRelationship # log or linear
mA = ds.XRayTubeCurrent
kVp = int(ds.KVP)

# plot the image
fig, ax = plt.subplots(1,1,figsize=[8,6])
ax.imshow(M, cmap='gray')
ax.set_title(f'"{pir}" data, {cols} x {rows} pixels, {mA} mA, {kVp} kVp')
plt.show()
```



This is a image of the anthropomorphic phantom we scanned. The “LOG” indicates the pixel values are not raw and have been through some post-processing. There are also header attributes for the image size, tube current, and kVp.

Let’s confirm that all the images look reasonable and check whether they use raw or log pixel values.

```
[3]: yd, xd = 3,7
fig,ax = plt.subplots(yd,xd,figsize=[xd*2, yd*2])

for i in range(len(L)):
    ds = pydicom.dcmread(imd+L[i])


    # plot image
    ax[i//xd, i%xd].imshow(ds.pixel_array, cmap='gray')

    # print pixel relationship
    print(ds.PixelIntensityRelationship)

for i in range(xd*yd):
    # turn off axes so images look niver
    ax[i//xd, i%xd].axis('off')
```

```
fig.tight_layout()
plt.show()
```

```
LOG
LOG
LIN
LIN
LOG
LOG
LOG
LOG
LOG
LOG
LOG
LIN
LIN
LIN
LIN
LIN
LIN
LIN
```



This looks like all of our images, but some of the data is linear and some is logged. It is important to be mindful of this for our computations.

We can print out the entire Dicom header to see all available information.

```
[4]: print(ds)
```

```
Dataset.file_meta -----
(0002, 0000) File Meta Information Group Length  UL: 168
```

(0002, 0001) File Meta Information Version	OB: b'\x00\x01'
(0002, 0002) Media Storage SOP Class UID	UI: Digital X-Ray Image Storage
- For Processing	
(0002, 0003) Media Storage SOP Instance UID	UI:
1.2.840.113619.2.203.4.2147483647.1643241378.878071	
(0002, 0010) Transfer Syntax UID	UI: Explicit VR Little Endian
(0002, 0012) Implementation Class UID	UI: 1.2.840.113619.6.203

(0008, 0005) Specific Character Set	CS: 'ISO_IR 100'
(0008, 0008) Image Type	CS: ['ORIGINAL', 'PRIMARY', '']
(0008, 0016) SOP Class UID	UI: Digital X-Ray Image Storage
- For Processing	
(0008, 0018) SOP Instance UID	UI:
1.2.840.113619.2.203.4.2147483647.1643241378.878071	
(0008, 0020) Study Date	DA: '20220126'
(0008, 0021) Series Date	DA: '20220126'
(0008, 0022) Acquisition Date	DA: '20220126'
(0008, 0023) Content Date	DA: '20220126'
(0008, 002a) Acquisition DateTime	DT: '20220126T175618.000000'
(0008, 0030) Study Time	TM: '161150.000'
(0008, 0031) Series Time	TM: '161408.000'
(0008, 0032) Acquisition Time	TM: '175618.000'
(0008, 0033) Content Time	TM: '175618.000'
(0008, 0050) Accession Number	SH: ''
(0008, 0060) Modality	CS: 'DX'
(0008, 0068) Presentation Intent Type	CS: 'FOR PROCESSING'
(0008, 0070) Manufacturer	LO: 'GE Healthcare'
(0008, 0080) Institution Name	LO: 'UNIV. OF CHICAGO'
(0008, 0081) Institution Address	ST: '5835 S. COTTAGE
GROVE\r\nMITCHELL TC114\r\nCHICAGO\r\nIL\r\n60637'	
(0008, 0090) Referring Physician's Name	PN: ''
(0008, 1010) Station Name	SH: 'MITTC114'
(0008, 1030) Study Description	LO: ''
(0008, 103e) Series Description	LO: 'Thumb'
(0008, 1040) Institutional Department Name	LO: 'RADIOLOGY'
(0008, 1050) Performing Physician's Name	PN: ''
(0008, 1070) Operators' Name	PN: ''
(0008, 1090) Manufacturer's Model Name	LO: 'Discovery XR656'
(0008, 2218) Anatomic Region Sequence 1 item(s) ----	
(0008, 0100) Code Value	SH: 'T-D8810'
(0008, 0102) Coding Scheme Designator	SH: 'SNM3'
(0008, 0104) Code Meaning	LO: 'Thumb'

(0010, 0010) Patient's Name	PN: 'Radiography Lab'
(0010, 0020) Patient ID	LO: 'MPQC01262021'
(0010, 0030) Patient's Birth Date	DA: ''
(0010, 0032) Patient's Birth Time	TM: ''
(0010, 0040) Patient's Sex	CS: 'O'

(0010, 1010) Patient's Age AS: ''
 (0011, 0010) Private Creator LO: 'GEMS_GDXE_FALCON_04'
 (0011, 1003) [Processed Series UID] UI:
 1.2.840.113619.2.203.4.2147483647.1643235248.198669
 (0011, 1004) [Acquisition Type] CS: 'SE'
 (0011, 1005) [Acquisition UID] UI:
 1.2.840.113619.2.203.4.2147483647.1643241375.209945
 (0011, 1006) [Image Dose] DS: '0.372645'
 (0011, 1009) [Non-Digital Exposures] SL: 20
 (0011, 1010) [Total Exposures] SL: 27
 (0011, 1012) [Patient Size String] LT: 'MEDIUM_ADULT'
 (0011, 1013) [SPS UID] UI:
 1.2.840.113619.2.203.4.2147483647.1643235110.550923
 (0011, 1015) [Detector ARC Gain] DS: '4400.0'
 (0011, 1016) [Processing Debug Info] LT: '1100,\n{1634 1559 414 1569
 414 469 1634 479 }(8
 elements)\r\n1200,\r\n1201,\r\n2100,\r\n2500,\r\n2501,\r\n2502,\r\n'
 (0011, 1017) [Override mode] CS: 'DISABLE'
 (0011, 1019) [Film Speed Selection] DS: '100.0'
 (0011, 1031) [Detected Field of View] IS: [1632, 1547, 419, 1557,
 419, 449, 1632, 459]
 (0011, 1034) [Compensated Detector Exposure] DS: '55.839879'
 (0011, 1036) [Median Anatomy Count Value] DS: '12027.0'
 (0011, 1037) [DEI lower and upper limit values] DS: [-2.000000, -1.000000,
 1.000000, 2.000000]
 (0011, 1042) [View IP] LO: 'antero-posterior'
 (0011, 1044) [Receptor Type] CS: 'TABLE'
 (0011, 1047) [Unknown] DS: '810.058228'
 (0011, 1059) Private tag data CS: 'NORMAL'
 (0011, 1060) Private tag data CS: 'HORIZONTAL'
 (0011, 1064) Private tag data CS: 'URP_DETECTOR'
 (0011, 1065) Private tag data LO: 'GlobalURPTable'
 (0011, 1066) Private tag data SH: 'No'
 (0011, 1067) Private tag data SH: 'Front'
 (0011, 1068) Private tag data LO: '40:f4:a0:0:c:4c'
 (0011, 1069) Private tag data DS: '-1.0'
 (0011, 106d) Private tag data DS: '218.5712'
 (0011, 1076) Private tag data CS: ['R', 'H']
 (0011, 1077) Private tag data CS: 'YES'
 (0011, 1078) Private tag data DS: '128.0'
 (0011, 1079) Private tag data DS: '14152.05874'
 (0011, 107a) Private tag data SS: 8
 (0011, 107b) Private tag data SS: 8
 (0011, 107c) Private tag data CS: 'YES'
 (0011, 107e) Private tag data DS: '180.0'
 (0011, 1080) Private tag data LT: '2'
 (0011, 1081) Private tag data LO: 'AP'
 (0011, 1082) Private tag data LO: 'GRID_FREQ_70_LP_CM'

(0018, 0015) Body Part Examined	CS: 'THUMB'
(0018, 0060) KVP	DS: '125.0'
(0018, 1020) Software Versions	LO:
'dm_Platform_release-V14_FW48-4'	
(0018, 1030) Protocol Name	LO: 'THUMB_RIGHT'
(0018, 1110) Distance Source to Detector	DS: '903.0'
(0018, 1111) Distance Source to Patient	DS: '829.0'
(0018, 1147) Field of View Shape	CS: 'RECTANGLE'
(0018, 1149) Field of View Dimension(s)	IS: [243, 222]
(0018, 1150) Exposure Time	IS: '10'
(0018, 1151) X-Ray Tube Current	IS: '100'
(0018, 1152) Exposure	IS: '1'
(0018, 1153) Exposure in uAs	IS: '1008'
(0018, 115e) Image and Fluoroscopy Area Dose Pro	DS: '0.812122'
(0018, 1160) Filter Type	SH: 'NONE'
(0018, 1164) Imager Pixel Spacing	DS: [0.183610, 0.183610]
(0018, 1166) Grid	CS: ['FIXED', 'FOCUSED']
(0018, 1190) Focal Spot(s)	DS: '0.6'
(0018, 1405) Relative X-Ray Exposure	IS: '372'
(0018, 1411) Exposure Index	DS: '5887.733905'
(0018, 1412) Target Exposure Index	DS: '814.00001'
(0018, 1413) Deviation Index	DS: '8.6'
(0018, 1508) Positioner Type	CS: ''
(0018, 1700) Collimator Shape	CS: 'POLYGONAL'
(0018, 1720) Vertices of the Polygonal Collimato	IS: [1621, 1546, 401, 1556,
401, 456, 1621, 466]	
(0018, 5101) View Position	CS: 'AP'
(0018, 6000) Sensitivity	DS: '215.383703'
(0018, 7000) Detector Conditions Nominal Flag	CS: 'YES'
(0018, 7001) Detector Temperature	DS: '27.700001'
(0018, 7004) Detector Type	CS: 'SCINTILLATOR'
(0018, 7005) Detector Configuration	CS: 'AREA'
(0018, 700a) Detector ID	SH: 'UA45725-2'
(0018, 701a) Detector Binning	DS: [1, 1]
(0018, 7020) Detector Element Physical Size	DS: [0.200000, 0.200000]
(0018, 7022) Detector Element Spacing	DS: [0.200000, 0.200000]
(0018, 7024) Detector Active Shape	CS: 'RECTANGLE'
(0018, 7026) Detector Active Dimension(s)	DS: [409.600000, 409.600000]
(0018, 7030) Field of View Origin	DS: [0.000000, 0.000000]
(0018, 7032) Field of View Rotation	DS: '0.0'
(0018, 7034) Field of View Horizontal Flip	CS: 'NO'
(0018, 704c) Grid Focal Distance	DS: '1000.0'
(0018, 7050) Filter Material	CS: 'COPPER'
(0018, 7052) Filter Thickness Minimum	DS: '0.0'
(0018, 7054) Filter Thickness Maximum	DS: '0.0'
(0018, 7060) Exposure Control Mode	CS: 'MANUAL'
(0018, 7062) Exposure Control Mode Description	LT: ''
(0020, 000d) Study Instance UID	UI:

```

1.2.840.113619.2.203.4.2147483647.1643235110.551096
(0020, 000e) Series Instance UID          UI:
1.2.840.113619.2.203.4.2147483647.1643235248.198658
(0020, 0010) Study ID                     SH: 'Def_Study_ID'
(0020, 0011) Series Number                 IS: '1629'
(0020, 0013) Instance Number               IS: '7'
(0020, 0020) Patient Orientation            CS: ['R', 'H']
(0020, 0062) Image Laterality              CS: 'R'
(0028, 0002) Samples per Pixel             US: 1
(0028, 0004) Photometric Interpretation    CS: 'MONOCHROME1'
(0028, 0010) Rows                          US: 2022
(0028, 0011) Columns                       US: 2022
(0028, 0100) Bits Allocated                US: 16
(0028, 0101) Bits Stored                   US: 14
(0028, 0102) High Bit                     US: 13
(0028, 0103) Pixel Representation          US: 0
(0028, 0300) Quality Control Image         CS: 'NO'
(0028, 0301) Burned In Annotation          CS: 'NO'
(0028, 1040) Pixel Intensity Relationship   CS: 'LIN'
(0028, 1041) Pixel Intensity Relationship Sign SS: 1
(0028, 1052) Rescale Intercept             DS: '0.0'
(0028, 1053) Rescale Slope                 DS: '1.0'
(0028, 1054) Rescale Type                  LO: 'US'
(0028, 2110) Lossy Image Compression       CS: '00'
(0032, 1060) Requested Procedure Description LO: ''
(0040, 0244) Performed Procedure Step Start Date DA: '20220126'
(0040, 0245) Performed Procedure Step Start Time TM: '161150.000'
(0040, 0253) Performed Procedure Step ID    SH: '164323511090911'
(0040, 0254) Performed Procedure Step Descriptio LO: 'Performed Desc'
(0040, 0302) Entrance Dose                 US: 0
(0040, 0310) Comments on Radiation Dose    ST: '125%'
(0040, 0555) Acquisition Context Sequence 0 item(s) ----
(0040, 8302) Entrance Dose in mGy          DS: '0.372645'
(0045, 0010) Private Creator               LO: 'GEMS_FALCON_03'
(0045, 1062) [User Window Center]          IS: '7140'
(0045, 1063) [User Window Width]          IS: '14024'
(0045, 1065) [Requested Detector Entrance Dose] IS: '14250'
(0045, 1069) [Collimator rotation]        IS: '0'
(0045, 1072) [Collimator Width]           IS: '218'
(0045, 1073) [Collimator Height]          IS: '244'
(0054, 0220) View Code Sequence 1 item(s) ----
    (0008, 0100) Code Value                SH: 'R-10206'
    (0008, 0102) Coding Scheme Designator   SH: 'SNM3'
    (0008, 0104) Code Meaning               LO: 'antero-posterior'
    -----
(0088, 0200) Icon Image Sequence 1 item(s) ----
    (0028, 0002) Samples per Pixel          US: 1
    (0028, 0004) Photometric Interpretation CS: 'MONOCHROME1'

```

(0028, 0010) Rows	US: 64
(0028, 0011) Columns	US: 64
(0028, 0100) Bits Allocated	US: 8
(0028, 0101) Bits Stored	US: 8
(0028, 0102) High Bit	US: 7
(0028, 0103) Pixel Representation	US: 0
(7fe0, 0010) Pixel Data	OB: Array of 4096 elements

(2050, 0020) Presentation LUT Shape	CS: 'INVERSE'
(7fdf, 0010) Private Creator	LO:
'GEMS_GDXE_ATHENAV2_INTERNAL_USE'	
(7fdf, 1010) [PPS Stream]	LT: 'PPSTREAM XRPP: 2.0 FILE /database/PixelData///1.2.840.113619.2.203.4.2147483647.1643235110.551096/1.2.840.113619.2.203.4.2147483647.1643235248.198658/1.2.840.113619.2.203.4.2147483647.1643241378.878071'
(7fdf, 1011) [Pixel Data References (temporary)]	LT:
'PPSTREAM XRPP: 2.0 FILE /database/PixelData//DefaultImageLocation-8366199105'	
(7fdf, 1025) [PPS Status]	CS: 'COMPLETED'
(7fe0, 0000) Group Length	UL: 8176980
(7fe0, 0010) Pixel Data	OW: Array of 8176968 elements

[]: